



SPECIFICATIONS

Items	Specifications
Stand	Dual light stable stand with Ball-bearing focusing controls
Viewing Body	Binocular tube, 45°, IPD 54-74mm
Eyepiece	Wide Field eyepiece WF 10x/20mm focusable paired with foldable eye guards.
Zoom ratio	4.4 : 1
Mechanical Stage	8x-35x
Condenser	Sub-stage abbe condenser NA 1.25 with aspheric lens. Iris diaphragm with snap-in blue filter. Rack and pinion movements on metal guide.
Focusing	Co-axial coarse and fine focusing on ball drive system for smooth movement with adjustable Tension Control knob.
Illumination	Built in LED ring light with +/-30 degree oblique incident illumination settings, Rheostat on ring light provides for intensity control, LED base light for transmitted illumination.
Electrical	Universal power supply (180V-240V, 50/60 Hz) UL, CE approved.
Packing	PACKED in Styrofoam box, with operation manual, allen wrench, dust cover, cleaning cloth, power cord and spare fuses.



Marketed by: Superfit Continental Pvt. Ltd.

5/6 Old Sitaram Building, First Floor, 204 Princess Street,
Mumbai 400002 INDIA Contact : +91-22-22018954 | 22054361
Email:- info@superfitindia.com, Website : www.superfitindia.com