



SPECIFICATIONS

Items	Specifications
Stand	Single mold Sturdy stand with anti rust materials. Extended base with hand rests for enhanced stability and comfort.
Viewing Body	Binocular or Trinocular tube 45° inclined, 360° rotatable, Interpupillary distance adjustable 54-74mm. Vacuum coated Prism for better transmission and anti fungus character.
Eyepiece	Focusable eyepiece with graduated scale WF10x/18mm F.O.V., with foldable eye guard, lockable, anti-fungus coating.
Nosepiece	Quadruple nosepiece (Ball bearing type) with rubber grip. PAR-CENTERED & PARFOCAL on CNC Machines.
Objectives	LP series DIN Plan Achromatic objectives 4x, 10x, 40x (spring loaded), 100x Oil (spring loaded) anti-fungus coating.
Mechanical Stage	Mechanical stage 135x125mm with co-axial drive, highly smooth on ball slides, Stainless Steel Holder for slide 75 x 50mm movement.
Condenser	Sub-stage abbe condenser NA 1.25 with aspheric lens. Iris diaphragm with snap-in blue filter. Rack and pinion movements on metal guide.
Focusing	Co-axial coarse and fine focusing on ball drive system for smooth Movement.
Illumination	Halogen 6V-20W illumination with variable illumination control.
	OR
Illumination	LED illumination with variable illumination control. Rechargeable Lithium Ion Battery Provide long Hours of Battery Backup incase of electricity failure.
Electrical	For HL Input 220v-240v AC & For LED Input 180v - 250v AC with removable imported cord.
Packing	PACKED in Styrofoam box, with operation manual, dust cover, cleaning cloth, power cord.
Optional Accessories	Phase Contrast attachment, Dark Field attachment(Wet Condenser), Simple Polarising attachment.



Bino



Trino



360°



Transmitted



Abbe



Polarized



PLAN



1 YEARS
Warranty



Phase



Dark



LED



Battery



Power Cord



Halogen



Marketed by: Superfit Continental Pvt. Ltd.

5/6 Old Sitaram Building, First Floor, 204 Princess Street,
Mumbai 400002 INDIA Contact : +91-22-22018954 | 22054361
Email:- info@superfitindia.com, Website : www.superfitindia.com

SUPERFIT®